

Follow this link to view the full presentation,
with animated (scans)

https://docs.google.com/presentation/d/10Z646nL0SNAAtWqknQ7BCpWcpEvJxT_GHyK5vZl9EGXo

Another **Philips Pulseq Interpreter**

from basics to pTx


ISMRM Virtual Meeting - 17 November 2023
Thomas Roos



UMC Utrecht

Pulseq is like music


While Philips
invented CDs...



Philips MRI
scanners
are *no* CD players

Why is Pulseq challenging?

Philips platform is *too* optimised - or is it?



Aim

Compatible


→ Current sequences & timing

Generic & easy



- No conversions
- Runs .seq directly
- No required labels, etc

Accurate

→ No timing changes



Architecture



Runtime



Modifies dummy sequence to convert Pulseq blocks to Philips 'blocks'




Images



Philips encoding

Recon




Raw




Pulseq encoding

Pulseseq - EPI		Voxel	Tra	Rel. SNR	TE	TR		
< 00:01		-	-	-	2.1	200		
Initial	Geometry	Contrast	Motion	Dyn/Ang	Postproc	Offc/Ang	Coils	Confl
+Pulseseq				yes				
+File [patch/#.seq]				epi_rs_noFS_ovs_phtime				
->Select (Edit it)								
1				1epi_rs_noFS_ovs_phtime				
2				6Nov3_2.5khz_SENS2x2_TE29...				
3				6nov_pTx1x1_TE296_TR90_Ts3...				
4				7nov_pTxONLY				

Playing: *a gradient echo*




Pulseq GRE
in Philips simulator




Pulseq & Native GRE
using Philips Achieva 7T
version R5.4


Playing: a spiral readout



Requested
in .seq




Measured
using Skope



k-space
trajectory

Playing: *an EPI ?!*




Playing: *a MP-RAGE*


Philips
Achieva 7T

1mm³
resolution


2x2 SENSE
3min




Philips
sequence



Pulseseq
sequence



Let's make Universal Pulses, **Truly Universal!**



2D GRE
Single transmit channels

pTx-Pulseq


Full pTx

- Magnitude & phase

Backwards compatible

- No changes to MATLAB toolbox
- Current interpreter reads ch1

Repeating time-points for the different channels






```
tx_pattern = [ 1 0 0 0 0 0 0 0 ]; % Only Tx #1  
rf_singleTx = rf; % Start with regular RF pulse  
rf_singleTx.signal = reshape(rf.signal' .* tx_pattern, 1, []);  
rf_singleTx.t = repmat(rf.t, 1, num_tx);
```


Hybrid sequence

Best of both worlds

Pulseq
Features
Flexibility

Native
Optimizations
Recon
Accessibility




Playing: *Hybrid pTx-Pulseq*

1mm³
resolution

R=4
using CS

On-scanner
reconstruction



Safety checks

During scan:


- SAR (average power)
- System specs (slew, strength)
- System safety (grad temp)

Using Pulseq:

- PNS (MATLAB model)
- Gradient resonances

To be implemented:


- SAR prediction
- Gradient temp prediction



Current state

Highlights

- Timing flexibility
- Full pTx
- Hybrid sequences



Next steps

- Improve simulations
- Add labels to data
- On-scanner recon of Pulseq encoding
- Off-center FOV

Acknowledgement



PHILIPS

Edwin Versteeg
Kyungmin Nam
Mark Gosselink
Hans Hoogduin
Dennis Klomp
Jeroen Siero
Jannie Wijnen

and the rest!

